

Mobile handwashing station

The CliniWash

Position near any mains electricity socket. The CliniWash needs no water supply.



Essential for use anywhere an effective hand wash is required: the CliniWash mobile handwash station is entirely self contained and needs no mains water service or drainage connections.

Highly mobile, robust and easy to maintain

Features of the CliniWash handwash station:-

- Highly mobile smooth running lockable castors enable speedy deployment.
- Robust but lightweight aluminium frame construction.
- · Hands free operation under warm running water.
- It is entirely self contained, with the TEALtainer water container holding 9 litres of fresh water.
- The waste water is contained within a heavy duty polyethytene bag - when the fresh is used up the waste bag will be full.
- The CliniWash mobile handwash station provides up to 40 wet/soap/rinse cycles per filling of the TEALtainer.
- The CliniWash has nothing to adjust. Both the water temperature (approx. 38 - 40°C) and the wet/soap/rinse cycles are factory set.
- Minimal lifting: the superglide castors let you just wheel the CliniWash mobile handwash station to the nearest tap for filling and the nearest sink or drain for emptying.
- · Simple daily hygiene routine.
- The CWP230's wash and rinse cycle is designed to match the World Health Organization's Guidelines Clean hands protect against infection.

The soap and towel holders are not supplied.

Specification

Overall height 1290mm
Height to bowl 890mm
Width 455mm
Weight (dry) 18.9kg
Power supply 230V mains

Construction

- Rigid, lightweight aluminium frame.
- Smooth running lockable castors.
- The bowl is white acrylic capped ABS.
- Comes complete with 13 amp moulded plug.
- TEALtainer securely retains waste water in a separate heavy duty LLDPE polythene 9 litre bag.

Included with the unit

TEALtainer water container.





Unit 2
Chelmsley Wood Industrial Estate
Waterloo Avenue
Chelmsley Wood
Birmingham, B37 6QQ, UK
Tel: +44 (0) 121 770 0593
Email: enquiries@tealwash.com

Order online at www.tealwash.com