



# Clean hands save lives!



## HMW230 Hygienius MediWash

**Matches World Health Organization Guidelines**  
*"Clean hands protect against infection"*

The automatic sensor starts a pre-programmed hand wash cycle when hands are placed within range.



The  
Hygienius  
MediWash can be  
plugged in to any  
convenient electrical  
socket

Suitable for medical, care home  
and healthcare settings



The Hygienius MediWash is entirely self contained and needs no mains water service or drainage connections.

# Pre-programmed hand wash cycle

The MediWash wash / rinse cycle is designed to match the World Health Organization Guidelines "Clean hands protect against infection".

To start the timed wash programme:-

1. Place the hands by the Greensensor and it will flash.
2. Water runs for 4 seconds to enable the wetting of hands.
3. The water stops for 20 seconds to allow for the application of liquid soap and performance of the 9 step hand wash routine.
4. Appropriate hand washing techniques are displayed on the screen.
5. Hot water runs for 10 seconds to enable the soap to be rinsed off the hands.
6. The total wash time is approximately 35 seconds.

## MediWash features:-

- Hands free operation under warm running water.
- The water always at the optimum temperature of 40° - 43°C.
- Highly portable - the Hygienius MediWash can be plugged in to any convenient electrical socket.
- Nothing to reset or adjust - simple daily hygiene routine.
- Highly efficient use of water - up to 80 wash cycles per filling.
- The unit comes complete with two easy to use TEALTainers which ensures that fresh and waste water are kept separate and hygienically contained.

## Specification

- Overall height 1370 mm
- Height to bowl 880 mm
- Width 620 mm
- Depth 590 mm
- Weight (dry) 30.8 kg (without dispensers)
- TEALTainer weight 1.5 kg (empty), 11.5 kg (full)
- Power supply Mains electricity

## Construction

- The body is fully moulded in blue fleck high impact polyethylene.
- The bowl is white acrylic capped ABS.
- Comes complete with moulded electrical plug.
- TEALTainers securely retain waste water in a separate heavy duty LLDPE polythene bags.

## Included with the Hygienius MediWash

- 2 TEALTainer water containers.
- The backboard designed to accept the customers' own choice of towel and soap dispenser.



CHEVAUNNESINCLAIR

"Consider it done"

22B Old Hope Road, Kingston 5, Jamaica W.I.

Find us on: WhatsApp at 1-876-440-4400.

Call us at: 1-876-440-4400 & 1-876-395-4436.

Email us at: dr.sinclair@live.com.

Fax us at: 1-876-906-7423.